

Neo G-418 Sulfate, Solution, 100 % Activity (50 mg/ml)

#Cat: NB-58-0123B

Size: 10ml

#Cat: NB-58-0123H

Size: 100ml

Product Information

Neo G-418 Sulfate Solution, 100 % Activity (50 mg/ml) Cat. No. NB-58-0123B (10 ml), NB-58-0123H (100 ml)

General Information

Neo G-418 is used in the selection and maintenance of eukaryotic cells, stable transfected with neomycin resistance genes. Neo G-418 is an aminoglycoside antibiotic, related to Gentamicin, and exhibits toxicity towards both eukaryotic and prokaryotic cells. It is produced by *Micromonospora rhodorangea* and acts by binding the ribosome, thus inhibiting protein synthesis in both prokaryotic and eukaryotic cells.

Product Specifications

Appearance	Clear frozen liquid
CAS No.	108321-42-2
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened, store at $+4^{\circ}\text{C}$ and use within 4-6 weeks
Shipping conditions	Frozen (Dry Ice)
Thawing	Overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize.
Working concentration	Recommended final concentration (0.1 – 1.0 mg/ml) depending on the cell type: <ul style="list-style-type: none">• HeLa: 200 – 600 $\mu\text{g}/\text{ml}$• 3T3 cells: 500 – 1000 $\mu\text{g}/\text{ml}$• CHO: 200 – 400 $\mu\text{g}/\text{ml}$• HEK 293: 500 – 800 $\mu\text{g}/\text{ml}$• Jurkat cells: 600 – 700 $\mu\text{g}/\text{ml}$

Instructions for Use

- Do not use Neo G-418 with antibiotic/antifungal preparations (e.g. Pen/Strep). These agents are competitive inhibitors of Neo G-418. Other antibiotics are potentially cross-reactive as well.
- Good laboratory practice requires optimal concentration of biologically active Neo G-418 to select and maintain cells. This must be determined for each set of growth conditions. Neo G-418 is used in the concentration range of 100 – 200 $\mu\text{g}/\text{ml}$ for bacteria, or 200 – 500 $\mu\text{g}/\text{ml}$ for most mammalian cells. Concentrations of Neo G-418 required for maintenance of selected cell lines are typically $\leq 50\%$ compared to selection.
- It is recommended, whenever experimental conditions are altered, the optimal concentration of the product should be re-evaluated.

Precautions and Disclaimer

This product is for research use only

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (www.neo-biotech.com) or phone (+33 9 77 40 09 09).