

# CANDOR<sup>®</sup>

The ELISA Experts



**Distributed by:**

**CliniSciences Group**



# ASSAY DEFENDER®

## Total HAMA and interference blocker

Ready-to use assay diluent for reduction of interference caused by HAMA (human anti-mouse antibodies), rheumatoid factors, nonspecific binding, cross-reactivities and matrix effects.

Interference in immunoassays can cause false positive and false negative results in routine testing. Reasons for the increasing prevalence of different kinds of interference are for example:

- population aging
- a growing number of patients with:
  - chronic autoimmune diseases
  - allergic and inflammatory diseases
- a better limit of detection in modern assays



Interferences can be classified due to their biochemical reasons and impact on assay performance.

- 1 Interference caused by antibodies from patient samples**  
e.g. HAMA (human anti-mouse antibodies), HAAA (human anti-animal antibodies), heterophilic antibodies and rheumatoid factors from patient samples
- 2 Interference caused by endogenous components of the sample**  
e.g. albumins, complement, lysozyme, fibrinogen,  $\alpha$ -1 Antitrypsin, atypically high lipid-, salt- or sugar concentrations as well as atypical viscosities
- 3 Interference caused by assay components**  
Assay components - like fluorescent or enzymatic labels - can cross-react with substances from the sample or change binding properties of the assay antibodies




Fig. A: Normal sandwich assay with no interference → correct result

Fig. B: An interfering antibody forms a bridge between capture and detection antibody → false positive result

Fig. C: Binding of the interfering antibodies to the capture antibody → false negative result

Fig. D: Components of Assay Defender® prevent binding of the interfering antibody → correct result

**Assay Defender® enables effective elimination of all kinds of interference and helps to prevent false positive or false negative results.**

**Distributed by:**


**CliniSciences Group**



## Use of Assay Defender®

Assay Defender® is particularly suitable for human or animal specimens containing HAMA (human anti-mouse antibodies) or other AAA (anti-animal antibodies).

Assay Defender® is a ready-to-use dilution buffer for human or animal body fluid specimens. The previously used dilution buffer can be easily replaced by Assay Defender®. Calibrators must also be diluted in Assay Defender®. Dilution in Assay Defender® can be 1:2 or higher. A useful dilution factor for most applications is 1:10 (1 part specimen + 9 parts Assay Defender®).



### **Assay Defender® excels at preventing false positive results in interference samples**

Serum samples from patients with rheumatoid illnesses (Sample 1 and 2) were analyzed using an established research ELISA for the detection of hAxI, a biomarker for hepatocellular carcinoma. Sample 3 is an interference-free sample from a healthy donor. Samples were diluted either in pure PBS-T/BSA, PBS-T/BSA supplemented with a commercially available HAMA Blocker (HAMA Blocker 3), commercially available ready-to-use HAMA blocking diluents (HAMA Blocker 1 and HAMA Blocker 2), or CANDOR's Assay Defender®. Assay Defender® alone can completely eliminate interferences and hence give a true negative result (< 25 ng/mL) in both interference samples. The results of interference free samples are not affected by Assay Defender®.

#### order numbers

Assay Defender®	50 mL	180 050
Assay Defender®	125 mL	180 125
Assay Defender®	500 mL	180 500

Bulk quantities available on request for kit manufacturers



# CliniSciences Group

## Austria

Company: CliniSciences GmbH  
Address: Sternwartestrasse 76, A-1180  
Wien - Austria  
Telephone: +43 720 115 580  
Fax: +43 720 115 577  
Email: [oesterreich@clinisciences.com](mailto:oesterreich@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Belgium

Company: CliniSciences S.R.L  
Address: Avenue Stalingrad 52, 1000  
Brussels - Belgium  
Telephone: +32 2 31 50 800  
Fax: +32 2 31 50 801  
Email: [belgium@clinisciences.com](mailto:belgium@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Denmark

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [danmark@clinisciences.com](mailto:danmark@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Finland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [suomi@clinisciences.com](mailto:suomi@clinisciences.com)  
Web: <https://www.clinisciences.com>



## France

Company: CliniSciences S.A.S  
Address: 74 Rue des Suisses, 92000  
Nanterre- France  
Telephone: +33 9 77 40 09 09  
Fax: +33 9 77 40 10 11  
Email: [info@clinisciences.com](mailto:info@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Germany

Company: Biotrend Chemikalien GmbH  
Address: Wilhelm-Mauser-Str. 41-43,  
50827 Köln - Germany  
Telephone: +49 221 9498 320  
Fax: +49 221 9498 325  
Email: [info@biotrend.com](mailto:info@biotrend.com)  
Web: <https://www.biotrend.com>



## Iceland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [island@clinisciences.com](mailto:island@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Ireland

Company: CliniSciences Limited  
Address: Ground Floor, 71 lower Baggot street  
Dublin D02 P593 - Ireland  
Telephone: +353 1 6971 146  
Fax: +353 1 6971 147  
Email: [ireland@clinisciences.com](mailto:ireland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Italy

Company: CliniSciences S.r.l  
Address: Via Maremmana inferiore 378  
Roma 00012 Guidonia Montecelio - Italy  
Telephone: +39 06 94 80 56 71  
Fax: +39 06 94 80 00 21  
Email: [italia@clinisciences.com](mailto:italia@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Netherlands

Company: CliniSciences B.V.  
Address: Gaetano Martinolaan 85, 6229  
GS Maastricht - Netherlands  
Telephone: +31 85 2082 351  
Fax: +31 85 2082 353  
Email: [nederland@clinisciences.com](mailto:nederland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Norway

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [norge@clinisciences.com](mailto:norge@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Poland

Company: CliniSciences sp.Z.o.o.  
Address: ul. Rotmistrza Witolda Pileckiego 67  
lok. 200 - 02-781 Warszawa -Poland  
Telephone: +48 22 307 0535  
Fax: +48 22 307 0532  
Email: [polska@clinisciences.com](mailto:polska@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Portugal

Company: Quimigen Unipessoal LDA  
Address: Rua Almada Negreiros, Lote 5, Loja 14,  
2615-275 Alverca Do Ribatejo - Portugal  
Telephone: +351 30 8808 050  
Fax: +351 30 8808 052  
Email: [info@quimigen.com](mailto:info@quimigen.com)  
Web: <https://www.quimigen.pt>



## Spain

Company: CliniSciences Lab Solutions  
Address: C/ Hermanos del Moral 13  
(Bajo E), 28019, Madrid - Spain  
Telephone: +34 91 269 40 65  
Fax: +34 91 269 40 74  
Email: [espana@clinisciences.com](mailto:espana@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Sweden

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [sverige@clinisciences.com](mailto:sverige@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Switzerland

Company: CliniSciences Limited  
Address: Marktgasse 18 8302 Kloten -  
Switzerland  
Telephone: +41 (044) 805 76 81  
Fax: +41 (044) 805 76 75  
Email: [switzerland@clinisciences.com](mailto:switzerland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## UK

Company: CliniSciences Limited  
Address: 11 Progress Business center, Whittle  
Parkway, SL1 6DQ Slough- United Kingdom  
Telephone: +44 (0)1753 866 511  
or +44 (0) 330 684 0982  
Fax: +44 (0)1753 208 899  
Email: [uk@clinisciences.com](mailto:uk@clinisciences.com)  
Web: <https://www.clinisciences.com>



## USA

Company: Biotrend Chemicals LLC  
Address: c/o Carr Riggs Ingram,  
500 Grand Boulevard, Suite 210 Miramar  
Beach, FL 32550- USA  
Telephone: +1 850 650 7790  
Fax: +1 850 650 4383  
Email: [info@biotrend-usa.com](mailto:info@biotrend-usa.com)  
Web: <https://www.biotrend-usa.com>

