

## FIX&PERM<sup>®</sup> Cell Fixation and Permeabilization

### Reaching targets for better flow cytometry staining

Intracellular staining of protein biomarkers using antibodies in flow cytometry now allows a number of ways for phenotypic and functional characterization of cells in both research and clinical diagnosis.

Definitive phenotypic identification of certain cell types can require labelling of both cell surface and intracellular markers. The cytoplasmic localization of typical membrane molecules can also prove to be a reliable lineage marker, such as cytoplasmic CD3 and CD22 in undifferentiated leukemia.

We can provide cell fixation and permeabilization reagents which assist antibodies to optimally stain their cell protein biomarkers.

### FIX&PERM<sup>®</sup>



- Mildly fixes cells, preserving their flow cytometric scatter characteristics
- Allows simultaneous characterization of both intracellular and cell surface markers
- Rapid technique – whole procedure can be carried out in less than one hour, ready for immediate analysis or storage for 24 hours at 4°C
- Stringent QC procedures – the quality of each lot is determined using well-defined blood samples and subsequent comparison of scatter characteristics of obtained leukocyte populations, ensuring consistent and reliable results lot after lot
- A range of intracellular antibodies with optimized protocols for use with FIX&PERM<sup>®</sup>

Choose from our range of complete kits or separate components:

Product Name	Quantity	Catalog Number
FIX&PERM <sup>®</sup> Kit (CE)	2x20ml	orb180078
FIX&PERM <sup>®</sup> Sample Kit (CE)	2x5ml	orb180079
FIX&PERM <sup>®</sup> Kit 1000 (CE)	2x100ml	orb180080
FIX&PERM <sup>®</sup> Solution A (Fix) (CE)	1x100ml	orb180081
FIX&PERM <sup>®</sup> Solution B (Perm) (CE)	1x100ml	orb180082
NM LYSE: Flow Cytometry Lysing	30ml	orb180088
NM LYSE: Flow Cytometry Lysing Solution	100ml	orb180089

### How does it work?

FIX&PERM<sup>®</sup> is a simple procedure making use of two reagents. Reagent A gently fixes cells, while Reagent B permeabilizes them. The specific formulations reduce background and allow simultaneous addition of permeabilization reagent and fluorochrome labelled antibodies.





## Products optimized for use with FIX&PERM® for leukemia detection via flow cytometry

The majority of available monoclonal antibody conjugates can be used with FIX&PERM® and are CE certified. However, in some cases the fixation step may have to be optimized.

We offer a number of mouse anti-human conjugated (PE/FITC) monoclonal antibodies ready for use with FIX&PERM®.

Antigen	Clone	Isotype	Format	Quantity	Catalog Number
CD3	UCHT1	IgG1	Purified	0.2 mg	orb180041
CD3	UCHT1	IgG1	PE-labeled	100 tests	orb180043
CD14	MEM18	IgG1	Purified	0.2 mg	orb180061
CD14	MEM18	IgG1	FITC-labeled	100 tests	orb180062
CD22	RFB4	IgG1	Purified	0.2 mg	orb180044
CD22	RFB4	IgG1	FITC-labeled	100 tests	orb180045
CD22	RFB4	IgG1	PE-labeled	100 tests	orb180046
CD45	VIT200	IgG2a	Purified	0.2 mg	orb180064
CD45	VIT200	IgG2a	FITC-labeled	100 tests	orb180065
CD65s	VIM2	IgM	Purified	0.2 mg	orb180066
CD65s	VIM2	IgM	FITC-labeled	100 tests	orb180067
CD68 (Macrosialin)	Ki-M7	IgG1	FITC-labeled	100 tests	orb180047
Lactoferrin	4C5	IgG1	Purified	0.2 mg	orb180051
Lactoferrin	4C5	IgG1	PE-labeled	100 tests	orb180052
Myeloperoxidase	8E2	IgG1	Purified	0.2 mg	orb180048
Myeloperoxidase	8E2	IgG1	FITC-labeled	100 tests	orb180049
Myeloperoxidase	8E2	IgG1	PE-labeled	100 tests	orb180050
Lysozyme	LZ-2	IgG1	Purified	0.2 mg	orb180053
Lysozyme	LZ-2	IgG1	FITC-labeled	100 tests	orb180054
Lysozyme	LZ-2	IgG1	PE-labeled	100 tests	orb180055
Cd68	Ki-M7	IgG1	FITC-labeled	100 tests	orb180047
Negative Control	VI-AP	IgG1	Purified	100 tests	orb180083
Negative Control	VI-AP	IgG1	FITC-labeled	100 tests	orb180084
Negative Control	VI-AP	IgG1	PE-labeled	100 tests	orb180085

Furthermore, a number of Combination Panels are offered for dual immunostaining protocols.

Antigen	Clone	Isotype	Format	Quantity	Catalog Number
Negative control	UCHT1	IgG1	FITC and PE	50 tests	orb180086
IgG2a Negative Control and IgG1 Negative Control	UCHT1	IgG1 and IgG2a	FITC and PE	50 tests	orb180087
Myeloperoxidase-C2 and Lactoferrin	UCHT1	IgG1	FITC and PE	50 tests	orb180056
Myeloperoxidase-C2 and CD3	MEM18	IgG1	FITC and PE	50 tests	orb180057
Myeloperoxidase-C2 and CD22	MEM18	IgG1	FITC and PE	50 tests	orb180058
Lysozyme and Lactoferrin	RFB4	IgG1	FITC and PE	50 tests	orb180059
CD45 and CD14	RFB4	IgG1 and IgG2a	FITC and PE	50 tests	orb180060

**Distributed by:**

**CliniSciences Group**

# CliniSciences Group

## Austria

Company: CliniSciences GmbH  
Address: Sternwartestrasse 76, A-1180  
Wien - Austria  
Telephone: +43 720 115 580  
Fax: +43 720 115 577  
Email: [oesterreich@clinisciences.com](mailto:oesterreich@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Belgium

Company: CliniSciences S.R.L  
Address: Avenue Stalingrad 52, 1000  
Brussels - Belgium  
Telephone: +32 2 31 50 800  
Fax: +32 2 31 50 801  
Email: [belgium@clinisciences.com](mailto:belgium@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Denmark

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [danmark@clinisciences.com](mailto:danmark@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Finland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [suomi@clinisciences.com](mailto:suomi@clinisciences.com)  
Web: <https://www.clinisciences.com>



## France

Company: CliniSciences S.A.S  
Address: 74 Rue des Suisses, 92000  
Nanterre- France  
Telephone: +33 9 77 40 09 09  
Fax: +33 9 77 40 10 11  
Email: [info@clinisciences.com](mailto:info@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Germany

Company: Biotrend Chemikalien GmbH  
Address: Wilhelm-Mauser-Str. 41-43,  
50827 Köln - Germany  
Telephone: +49 221 9498 320  
Fax: +49 221 9498 325  
Email: [info@biotrend.com](mailto:info@biotrend.com)  
Web: <https://www.biotrend.com>



## Iceland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [island@clinisciences.com](mailto:island@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Ireland

Company: CliniSciences Limited  
Address: Ground Floor, 71 lower Baggot street  
Dublin D02 P593 - Ireland  
Telephone: +353 1 6971 146  
Fax: +353 1 6971 147  
Email: [ireland@clinisciences.com](mailto:ireland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Italy

Company: CliniSciences S.r.l  
Address: Via Maremmana inferiore 378  
Roma 00012 Guidonia Montecelio - Italy  
Telephone: +39 06 94 80 56 71  
Fax: +39 06 94 80 00 21  
Email: [italia@clinisciences.com](mailto:italia@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Netherlands

Company: CliniSciences B.V.  
Address: Kraaijenhoffstraat 137A,  
1018RG Amsterdam, Netherlands  
Telephone: +31 85 2082 351  
Fax: +31 85 2082 353  
Email: [nederland@clinisciences.com](mailto:nederland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Norway

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [norge@clinisciences.com](mailto:norge@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Poland

Company: CliniSciences sp.Z.o.o.  
Address: ul. Rotmistrza Witolda Pileckiego 67  
lok. 200 - 02-781 Warszawa -Poland  
Telephone: +48 22 307 0535  
Fax: +48 22 307 0532  
Email: [polska@clinisciences.com](mailto:polska@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Portugal

Company: Quimigen Unipessoal LDA  
Address: Rua Almada Negreiros, Lote 5, Loja 14,  
2615-275 Alverca Do Ribatejo - Portugal  
Telephone: +351 30 8808 050  
Fax: +351 30 8808 052  
Email: [info@quimigen.com](mailto:info@quimigen.com)  
Web: <https://www.quimigen.pt>



## Spain

Company: CliniSciences Lab Solutions  
Address: C/ Hermanos del Moral 13  
(Bajo E), 28019, Madrid - Spain  
Telephone: +34 91 269 40 65  
Fax: +34 91 269 40 74  
Email: [espana@clinisciences.com](mailto:espana@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Sweden

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [sverige@clinisciences.com](mailto:sverige@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Switzerland

Company: CliniSciences Limited  
Address: Marktgasse 18 8302 Kloten -  
Switzerland  
Telephone: +41 (044) 805 76 81  
Fax: +41 (044) 805 76 75  
Email: [switzerland@clinisciences.com](mailto:switzerland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## UK

Company: CliniSciences Limited  
Address: 11 Progress Business center, Whittle  
Parkway, SL1 6DQ Slough- United Kingdom  
Telephone: +44 (0)1753 866 511  
or +44 (0) 330 684 0982  
Fax: +44 (0)1753 208 899  
Email: [uk@clinisciences.com](mailto:uk@clinisciences.com)  
Web: <https://www.clinisciences.com>



## USA

Company: Biotrend Chemicals LLC  
Address: c/o Carr Riggs Ingram,  
500 Grand Boulevard, Suite 210 Miramar  
Beach, FL 32550- USA  
Telephone: +1 850 650 7790  
Fax: +1 850 650 4383  
Email: [info@biotrend-usa.com](mailto:info@biotrend-usa.com)  
Web: <https://www.biotrend-usa.com>

