

---

**Mouse Anti -Human IgG2 antibody {FITC}**

#Cat: NB-22-31353-100 Size: 100µl

**General Information**

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	FITC-conjugated mouse monoclonal anti-Human IgG2 secondary antibody. For use in most research applications.
<b>Applications</b>	WB, IHC-P, IF, ICC, ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**Product Properties**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	4-5 mg/ml
<b>Conjugation</b>	FITC
<b>Purification</b>	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	WB 1:5000-20000 ELISA 1:10000-20000 IHC 1:5000-10000 IF 1:5000-10000
<b>Formulation</b>	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide
<b>Isotype</b>	IgG2
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**Target Information**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.