

---

## 5-FAM antibody

### Product Information

Catalog No.:	FNab10903
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.06% PC300 and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	KLH conjugated 5-FAM
Synonyms:	5-Carboxyfluorescein
Observed MW:	None
Uniprot ID :	None

### Application

Reactivity:	FAM
Tested Application:	ELISA
Recommended dilution:	1:5000-1:10000