

Mini Neo-Dry Bath Incubator N70 #Cat: NB-12-0011

The N70 Neo-Dry block heating is designed for testing milk for antibacterial substances such as antibiotics and sulphonamideresidues. It includes a digital thermometer and timer. Heater is either AC or DC compatible. AC adapter is included, DC adapter for vehicle must be purchased separately.

Features:

- A personal dry block heating and for digital dry block heating
- It pretty much fits anywhere on your lab bench
- It has a very simple touch pad digital display for temperature and time selection
- It's only 4-inch wide
- High quality -- low price



Specifications:

Model	N70
Cat#	NB-12-0011
Power supply	DC12V
Power	15 W
Average power consumption	5 W
Temperature range	RT+5°C ~ 80°C
Timing range	1h ~ 999h 1min ~ 999min 1sec ~ 999 sec (optional)
Accuracy of the temperature	≤± 0.5°C
Display accuracy	0.1°C
Temp. Homogeneity	±0.3°C
Temp. Stability	±0.3°C
Heating time (from 25°C to 80°C)	≤20min (The ambient temperature is 20°C~30°C)
Ambient temperature	5°C ~ 35°C
Dimensions(W×D×H)	99×135×34mm
Weight	0.5 kg