Product Description:	NT1 antibody, human, C-term. NTR1 Neurotensin 1 receptor G Protein-coupled receptor Polyclonal antibody raised against the carboxy human NT1 receptor.	y-terminus of the
Product-No.:	S12-9036-A-G	
Immunogen:	residues 398-418 coupled to KLH (NP_002522) KLH- kadsvssnhtlssnatretly	
Host:	Rabbit	
Form:	IgG fraction supplied in PBS pH 7.4 containing 0.05% thimerosal as preservative Antibodies were purified by peptide affinity chromatography	
Quantity:	50 µg	anti-NT1
Concentration:	0.1 mg/ml	EZ EZ 200 -
Applications:	IC, W W = Western blot (1 µg/ml) IC = Immunocytochemistry (1 µg/ml)	120 - 100 - 80 -
Species cross-reactivity:	human	70 - 60 -
Sensitivity:	50 ng	50 -
Storage:	4 °C short term -20 °C long term	40 -