

Monoclonal Anti-Cre Recombinase

#Cat: NB-01-0016 Size: 100 µl

Format: Crude Ascites fluid

Clones: 2CRE-2D8

Immunogen: Recombinant Cre recombinase (aa 77 to 343) as a fusion protein.

Murine Isotype: IgG1, k

Specificity: Cre recombinase. Epitope located between aa 77 and 343.

Applications:

ELISA

Western Blot (natural and recombinant proteins)

Immunofluorescence

Immunohistochemistry

Recommended dilutions: 1/500-1/5000

Storage:

Store at -20°C (up to 3 years)

Avoid repeat freezing and thawing cycles. Does not contain any preservative

These products are intended for research purposes only.

MATERIAL SAFETY DATA SHEET INFORMATION

To the best of our knowledge, these products do not require a Material Safety Data Sheet. However, all the properties of these products (and, if applicable, each of their components) have not been thoroughly investigated. Therefore, we recommend that you use gloves and eye protection, and wear a laboratory coat when working with these products.