

Anti-UGGT1 antibody (1456-1555) Cat # NB-22-20758

Size : 50 µL ; 100 µL ; 200 µL

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-UGGT1 (1456-1555) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Rat, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56886
Gene Symbol	UGGT1
Uniprot ID	UGGG1_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1456-1555 of human UGGT1 (NP_064505.1).
Immunogen Region	1456-1555
Specificity	
Immunogen sequence	



Western blot analysis of extracts of various cell lines, using UGGT1 antibody (STJ26034) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.