

Product datasheet PS118F

MONOSAN[®]

Polyclonal Goat anti-Rabbit IgG, IgM, IgA - FITC

Clone no. -

MONOSAN

Product name	Polyclonal Goat anti-Rabbit IgG, IgM, IgA - FITC
Host	Goat
Applications	WB, IHC (1:200), ELISA, IF, FACS, EM
Species reactivity	rabbit
Conjugate	FITC
Immunogen	IgG, IgM, IgA
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	1.2 mg/ml
Format	-
Storage buffer	0.05 M sodium phosphate, 0.1 M NaCl, 50% glycerol, 5 mg/ml BSA
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS118F

MONOSAN[®]

Polyclonal Goat anti-Rabbit IgG, IgM, IgA - FITC

Clone no. -

MONOSAN

Additional info

Solution 1.2 mg/mL antibody, 50% glycerol, 0.1M NaCl, 0.05M sodium phosphate pH 7.5, 5 mg/ml BSA. Rabbit immunoglobulins IgG, IgM, IgA

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES
