

Normal Bovine serum blocking buffer

#Cat: NB-23-00180-1

Size: 100mL

Intended use:

Storage: 2-8C

Normal Bovine Serum Blocking Buffer is used in immunohistochemical staining procedures to block non-specific staining when the host species of the secondary antibody is bovine (cattle).

Component:

| Catalog No | Size |
|---------------|----------------------|
| NB-23-00180-1 | 100mL (Ready to use) |

Solution contains 0.05% NaN₃.

Recommended protocol:

1. Apply 100uL or enough volume of **Normal Bovine Serum Blocking Buffer** on tissue before adding primary antibody and incubate for 10 minutes at room temperature.
2. Blot off Normal Bovine Serum Blocking Buffer Solution (Do not rinse) and add primary antibody.

Precaution:

Wear gloves and take other necessary laboratory safety procedure.

*For reference only For
Research Use Only.*