

A18721

Leader in Biomolecular Solutions for Life Science



C5b-9 Mouse mAb

Catalog No.: A18721

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

SMab Recombinant Monoclonal
Antibody

Applications

IHC-P,IF/ICC,ELISA

Cross-Reactivity

Human,Horse,Pig,Monkey,Baboon

Background

Recommended Dilutions

IHC-P	1:50 - 1:200
IF/ICC	1:50 - 1:200

Immunogen Information

Gene ID

727/729/730/731/735

Swiss Prot

Immunogen

Recombinant protein of Human C5b-9.

Synonyms

Contact

 www.abclonal.com

Product Information

Source

Mouse

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.