

Proteus Mini Purification Spin Column Pack

#Cat: NB-45-00062

Size: 25pc

Neo Biotech's Proteus mini purification spin columns are designed for small scale protein purifications such as those required for expression trials, solubility determination tests, screening, titrating and scouting studies. The columns fit all standard microfuges and allow you to process multiple samples in parallel.



Specification:

Sinter type:	Ultra High Density Polyethylene
Plastic construction:	Polypropylene
Pore size:	10-40 μm (mean pore size, 18 μm)
Maximum sample volume:	0.65 ml
Hold-up volume	< 5 μl
Maximum g force:	14,000 g
Typical spin times:	1 min for 0.65 ml sample at 14,000 g
Storage:	Store at RT (non-sterile)
Shelf-life:	24 months

Ordering Information:

Product	Units	Order Code
Proteus Mini Purification Spin Column Pack	25	NB-45-00062