

LTBR antibody

Product Information

Catalog No.: FNab04887

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months(Avoid repeated freeze / thaw cycles.)

Background

Receptor for the heterotrimeric lymphotoxin containing LTA and LTB, and for TNFS14/LIGHT. Promotes apoptosis via TRAF3 and TRAF5. May play a role in the development of lymphoid organs.

Immunogen information

Immunogen: lymphotoxin beta receptor(TNFR superfamily, member 3)

Synonyms: Tumor necrosis factor receptor superfamily member 3|Lymphotoxin-beta

receptor|Tumor necrosis factor C receptor|Tumor necrosis factor

receptor 2-related protein|Tumor necrosis factor receptor type III (TNF-

RIII, TNFR-III)|LTBR|D12S370|TNFCR|TNFR3|TNFRSF3

Observed MW: 47 kDa Uniprot ID: P36941

Application


Reactivity: Human, Mouse

Tested Application: ELISA, WB, IHC, FC


Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200

Image:





Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab04887(LTBR antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab04887(LTBR antibody) at dilution of 1:1000