

Pepsin

Catalog # 6011

For Research Use Only - Not Human or Therapeutic Use

| | |
|--------------|--|
| DESCRIPTION: | Pepsin |
| APPLICATION: | Use to extract collagens from tissue samples in acidic buffer. Keep collagen samples at 4°C to avoid denaturation. |
| QUANTITY: | 1 mg |
| FORM: | Lyophilized powder |
| SOURCE: | Porcine gastric mucosa |
| PURITY: | 2500-3500 units/mg protein (E1%/280) |
| STORAGE: | -20°C |
| STABILITY: | 1 year |
| NOTES: | N/A |
| REFERENCES: | N/A |