

COG6 antibody

Product Information

Catalog No.:	FNab01830
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for normal Golgi function.

Immunogen information

Immunogen:	component of oligomeric golgi complex 6
Synonyms:	Conserved oligomeric Golgi complex subunit 6 (COG complex subunit 6) Component of oligomeric Golgi complex 6 COG6 KIAA1134
Observed MW:	73 kDa
Uniprot ID :	Q9Y2V7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab01830(COG6 antibody) at dilution of 1:300