

SFTPA1 antibody

Product Information

Catalog No.:	FNab07796
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	surfactant protein A1
Synonyms:	Pulmonary surfactant-associated protein A1 (PSP-A, PSPA, SP-A, SP-A1) 35 kDa pulmonary surfactant-associated protein Alveolar proteinosis protein Collectin-4 SFTPA1 COLEC4 PSAP SFTP1 SFTPA SFTPA1B
Observed MW:	None
Uniprot ID :	Q8IWL2

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:100 - 1:1000
Image:	



Immunohistochemistry of paraffin-embedded human lung tissue slide using FNab07796 (SFTP1 antibody) at dilution of 1:800