

Goat Anti-Human IgM

Product Information

Catalog No.:	FNab09908
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Goat
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen: human IgM

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution	:WB: 1:1000-1:10000

Note: This product is for research use only, not for diagnostic or therapeutic use.