

RP10003LQ

Leader in Biomolecular Solutions for Life Science



www.abclonal.com

Recombinant Human STAM-binding protein/AMSH protein

Catalog No.: RP10003LQ **Recombinant**

Sequence Information

Species	Gene ID	Swiss Prot
<I>E. coli</I>	10617	O95630

Tags

His Tag

Synonyms

STAMBp;AMSH;MICCAP

Product Information

Source

<I>E. coli</I>

Purification

~ 90% by SDS-PAGE

Endotoxin

Formulation

20 mM Tris, 150 mM NaCl, 2 mM β ME, 10% Glycerol

Reconstitution

Background

AMSH belongs to the JAMM family of deubiquitinases that are Zn²⁺-dependent metalloproteases. AMSH is an endosome-associating DUB that regulates endosomal sorting by deubiquitination. AMSH specifically cleaves lysine 63-polyubiquitin chains.

Basic Information

Description

Bio-Activity

Storage


Store at -70°C. This product is stable at $\leq -70^\circ\text{C}$ for up to 1 year from the date of receipt. For optimal storage, aliquot into smaller quantities after centrifugation and store at recommended temperature. Avoid repeated freeze-thaw cycles. Avoid repeated freeze/thaw cycles.

Contact



www.abclonal.com

Validation Data



Recombinant human STAM-binding protein/AMSH protein was determined by SDS-PAGE with Coomassie Blue.