



Normal Equine (horse)
serum blocking buffer

NB-23-00184-1

Normal Equine (horse) serum blocking buffer

#Cat: NB-23-00184-1

Size: 100mL

Intended use:

Storage: 2- 8°C

Normal Equine (Horse) Serum Blocking Buffer is used in immunohistochemical staining procedures to block nonspecific staining when the host species of the secondary antibody is Equine (Horse).

Kit components:

Catalog Number	Size
NB-23-00184-1	100 ml (ready to use)

Solution contains 0.05% NaN₃.

Recommended protocol:

1. Apply 100uL or enough volume of Normal Equine (Horse) Serum Blocking Buffer on tissue before adding primary antibody and incubate for 10 minutes at room temperature.
2. Blot off Normal Equine (Horse) Serum Blocking Buffer Solution (Do not rinse) and add primary antibody.

Precaution:

Wear gloves and take other necessary laboratory safety procedure.

Storage:

Store at 2-8°C.

FOR RESEARCH USE ONLY

Neo-Biotech

74 rue des Suisses – 92000 Nanterre

05/07/2024