

---

## ALB antibody

### Product Information


Catalog No.:	FNab00969
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6F11
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	BSA
Synonyms:	Albumin BSA ALB
Observed MW:	67 kDa
Uniprot ID :	P02769

### Application

Reactivity:	Bovine
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:5000-1:30000
Image:	



Bovine serum tissue were subjected to SDS PAGE followed by western blot with FNab00969( BSA Antibody) at dilution of 1:4000