Ne Biotech

| Product: | Monoclonal Antibody to Apolipoprotein E4 | | | | | |
|---|--|--------------|---|------------------------|--------------------------------|--|
| Immunogen: | na | | | | | |
| Pack sizes: Cat. No. | 100 µg: NB-02-1464 | Host Animal: | | Mouse | Lot No. na | |
| | 1 mg: NB-02-1762 | Format: | | Protein A Purification | | |
| IgG Concentration: na | | | Buffer: pH7.2 PBS | | | |
| Working Concentration: na | | | Preservative: 0.09% Sodium Azide | | | |
| Storage: | rage: + 2 - 8 | | | | | |
| Recommended Applications: | | | | | | |
| The antibody is suitable for the development of immunoassays or immunoaffinity purification columns | | | | | | |
| Sensitivity: | | | Specifity: | | | |
| Apolipoprotein E4 (~3ng/ml) | | | The antibody exhibited crossreactivity for the following compounds: Apolipoprotein E2 (<1%) Apolipoprotein E3 (<1%) | | | |
| Related Products: | | | | | | |
| na, na, | | | | | | |
| This product as supplied is intended for research applications only, not for use in therapeutic or diagnostic applications without the expressed written authorization of Randox Life Sciences. | | | | | | |
| | | | | Information corr | rect at time of going to print | |