

Anti-HLA-A, B, C antibody [W6/32] Cat # NB-22-62345-25

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Mouse monoclonal antibody anti-HLA-A, B, C is suitable for use in Flow Cytometry, Western

Blot, Immunohistochemistry and Immunoprecipitation research applications.

Applications FC, WB, IHC, IP

Host/Source Mouse **Reactivity** Human

PRODUCT PROPERTIES

Clonality Monoclonal
Clone ID W6/32
Concentration 0.5 mg/mL
Conjugation Unconjugated

Purification

Dilution Range

Formulation PBS with 0.05% Proclin300

Isotype IgG2ak

Storage Instruction Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID 3106 Gene Symbol HLA-B

Uniprot ID HLAB HUMAN

Immunogen

Immunogen Region

Specificity

Immunogen

sequence

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.