

## 1. PRODUCT AND COMPANY IDENTIFICATION

### 1.1 Product identifier

Product name : Monomethyl auristatin E  
Catalog No. : HY-15162  
CAS No. : 474645-27-7

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Laboratory chemicals, manufacture of substances.

### 1.3 Details of the supplier of the safety data sheet

Company: MedChemExpress Europe  
Tel: +4686500910  
E-mail: eu.sales@medchemexpress.com

### 1.4 Emergency telephone number

European (EU) emergency number: 112  
American (USA) emergency number: 911  
CHEMTREC emergency number: North America +1-800-424-9300 (toll free)  
Asia emergency number: 110  
Australia emergency number: 000  
All numbers are available 24/7/365

## 2. HAZARDS IDENTIFICATION

### 2.1 Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008 [GHS/CLP]

Acute toxicity, oral and acute toxicity dermal (Category 1, 2),H300+H330

Germ cell mutagenicity (Category 1A, 1B),H340

Reproductive toxicity (Category 1A, 1B),H360

Specific target organ toxicity, single exposure (Category 1),H370

Specific target organ toxicity, repeated exposure (Category 1),H372

Hazardous to the aquatic environment, long-term hazard (Category 3),H412

### 2.2 GHS Label elements, including precautionary statements

Labeling according to Regulation (EC) No 1272/2008 [CLP]

Pictogram



Signal word Danger

Hazard statement(s)

H300+H330 Fatal if swallowed or if inhaled

H340 May cause genetic defects

H360 May damage fertility or the unborn child

H370 Causes damage to organs

---

H372 Causes damage to organs through prolonged or repeated exposure

H412 Harmful to aquatic life with long lasting effects

Precautionary statement(s)

Prevention

P260 Do not breathe dust/fume/gas/mist/vapours/spray.

P264 Wash hands thoroughly after handling

P270 Do not eat, drink or smoke when using this product.

P271 Use only outdoors or in a well-ventilated area.

P273 Avoid release to the environment.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P284 Wear respiratory protection.

Response

P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P320 Specific treatment is urgent (see on this label).

P321 Specific treatment (see on this label).

P330 Rinse mouth.

Storage

P405 Store locked up.

Disposal

P501 Dispose of contents/container to in accordance with local regulation.

### 3. COMPOSITION/INFORMATION ON INGREDIENTS

#### 3.1 Substances

Synonyms: MMAE; SGD-1010; Vedotin

Formula:  $C_{39}H_{67}N_5O_7$

Molecular Weight: 717.98

CAS No. : 474645-27-7

### 4. FIRST AID MEASURES

#### 4.1 Description of first aid measures

##### Eye contact

Remove any contact lenses, locate eye-wash station, and flush eyes immediately with large amounts of water. Separate eyelids with fingers to ensure adequate flushing. Promptly call a physician.

##### Skin contact

Rinse skin thoroughly with large amounts of water. Remove contaminated clothing and shoes and call a physician.

##### Inhalation

Immediately relocate self or casualty to fresh air. If breathing is difficult, give cardiopulmonary resuscitation (CPR). Avoid mouth-to-mouth resuscitation.

##### Ingestion

---

Wash out mouth with water; Do NOT induce vomiting; call a physician.

#### 4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2).

#### 4.3 Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

### 5. FIRE FIGHTING MEASURES

#### 5.1 Extinguishing media

##### Suitable extinguishing media

Use water spray, dry chemical, foam, and carbon dioxide fire extinguisher.

#### 5.2 Special hazards arising from the substance or mixture

During combustion, may emit irritant fumes.

#### 5.3 Advice for firefighters

Wear self-contained breathing apparatus and protective clothing.

### 6. ACCIDENTAL RELEASE MEASURES

#### 6.1 Personal precautions, protective equipment and emergency procedures

Use full personal protective equipment. Avoid breathing vapors, mist, dust or gas. Ensure adequate ventilation. Evacuate personnel to safe areas.

Refer to protective measures listed in sections 8.

#### 6.2 Environmental precautions

Try to prevent further leakage or spillage. Keep the product away from drains or water courses.

#### 6.3 Methods and materials for containment and cleaning up

Absorb solutions with finely-powdered liquid-binding material (diatomite, universal binders); Decontaminate surfaces and equipment by scrubbing with alcohol; Dispose of contaminated material according to Section 13.

### 7. HANDLING AND STORAGE

#### 7.1 Precautions for safe handling

Avoid inhalation, contact with eyes and skin. Avoid dust and aerosol formation. Use only in areas with appropriate exhaust ventilation.

#### 7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly sealed in cool, well-ventilated area. Keep away from direct sunlight and sources of ignition.

Recommended storage temperature: -20°C, sealed storage, away from moisture and light

\* The compound is unstable in solutions, freshly prepared is recommended.

Shipping at room temperature if less than 2 weeks.

#### 7.3 Specific end use(s)

No data available.

### 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

#### 8.1 Control parameters

##### Components with workplace control parameters

---

This product contains no substances with occupational exposure limit values.

## 8.2 Exposure controls

### Engineering controls

Ensure adequate ventilation. Provide accessible safety shower and eye wash station.

### Personal protective equipment

#### Eye protection

Safety goggles with side-shields.

#### Hand protection

Protective gloves.

#### Skin and body protection

Impervious clothing.

#### Respiratory protection

Suitable respirator.

#### Environmental exposure controls

Keep the product away from drains, water courses or the soil. Clean spillages in a safe way as soon as possible.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

### 9.1 Information on basic physical and chemical properties

Appearance	Solid
Odor	No data available
Odor threshold	No data available
pH	No data available
Melting/freezing point	>90 °C
Boiling point/range	873.5±65.0°C at 760 mmHg
Flash point	482.1±34.3°C
Evaporation rate	No data available
Flammability (solid, gas)	No data available
Upper/lower flammability or explosive limits	No data available
Vapor pressure	No data available
Vapor density	No data available
Relative density	No data available
Water Solubility	No data available
Partition coefficient	No data available
Auto-ignition temperature	No data available
Decomposition temperature	No data available
Viscosity	No data available
Explosive properties	No data available
Oxidizing properties	No data available

### 9.2 Other safety information

No data available.

## 10. STABILITY AND REACTIVITY

### 10.1 Reactivity

No data available.

### 10.2 Chemical stability

---

Stable under recommended storage conditions.

### 10.3 Possibility of hazardous reactions

No data available.

### 10.4 Conditions to avoid

No data available.

### 10.5 Incompatible materials

Strong acids/alkalis, strong oxidising/reducing agents.

### 10.6 Hazardous decomposition products

Under fire conditions, may decompose and emit toxic fumes.

Other decomposition products - no data available.

## 11. TOXICOLOGICAL INFORMATION

### 11.1 Information on toxicological effects

#### Acute toxicity

Classified based on available data. For more details, see section 2

#### Skin corrosion/irritation

Classified based on available data. For more details, see section 2

#### Serious eye damage/irritation

Classified based on available data. For more details, see section 2

#### Respiratory or skin sensitization

Classified based on available data. For more details, see section 2

#### Germ cell mutagenicity

Classified based on available data. For more details, see section 2

#### Carcinogenicity

IARC: No component of this product present at a level equal to or greater than 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at a level equal to or greater than 0.1% is identified as a potential or confirmed carcinogen by ACGIH.

NTP: No component of this product present at a level equal to or greater than 0.1% is identified as a anticipated or confirmed carcinogen by NTP.

OSHA: No component of this product present at a level equal to or greater than 0.1% is identified as a potential or confirmed carcinogen by OSHA.

#### Reproductive toxicity

Classified based on available data. For more details, see section 2

#### Specific target organ toxicity - single exposure

Classified based on available data. For more details, see section 2

#### Specific target organ toxicity - repeated exposure

Classified based on available data. For more details, see section 2

#### Aspiration hazard

Classified based on available data. For more details, see section 2

#### Additional information

---

This information is based on our current knowledge. However the chemical, physical, and toxicological properties have not been completely investigated.

## 12. ECOLOGICAL INFORMATION

### 12.1 Toxicity

No data available.

### 12.2 Persistence and degradability

No data available.

### 12.3 Bioaccumulative potential

No data available.

### 12.4 Mobility in soil

No data available.

### 12.5 Results of PBT and vPvB assessment

PBT/vPvB assessment unavailable as chemical safety assessment not required or not conducted.

### 12.6 Other adverse effects

No data available.

## 13. DISPOSAL CONSIDERATIONS

### 13.1 Waste treatment methods

#### Product

Dispose substance in accordance with prevailing country, federal, state and local regulations.

#### Contaminated packaging

Conduct recycling or disposal in accordance with prevailing country, federal, state and local regulations.

## 14. TRANSPORT INFORMATION

### 14.1 ADR/RID

UN proper shipping name : TOXIC SOLID,ORGANIC, N. O. S.

UN number : 2811

Class : 6.1

Packing group : I

Environmental hazards :no

### 14.2 IMDG

UN proper shipping name : TOXIC SOLID,ORGANIC, N. O. S.

UN number : 2811

Class : 6.1

Packing group : I

Environmental hazards :no

### 14.3 IATA

UN proper shipping name : TOXIC SOLID,ORGANIC, N. O. S.

UN number : 2811

---

Class : 6.1

Packing group : I

Environmental hazards :no

## 15. REGULATORY INFORMATION

This safety datasheet complies with the requirements of Regulation (EC) No. 453/2010.

### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

No data available.

### 15.2 Chemical safety assessment

Chemical safety assessment was not carried out for this product.

## 16. OTHER INFORMATION

Copyright 2023 MedChemExpress. The above information is correct to the best of our present knowledge but does not purport to be all inclusive and should be used only as a guide. The product is for research use only and for experienced personnel. It must only be handled by suitably qualified experienced scientists in appropriately equipped and authorized facilities. The burden of safe use of this material rests entirely with the user. MedChemExpress disclaims all liability for any damage resulting from handling or from contact with this product.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: [tech@MedChemExpress.com](mailto:tech@MedChemExpress.com)

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA