

Anti-Ets1 (C-terminal region)

#Cat: NB-01-0033 Size: 100 μl

Immunogen: synthetic peptide coupled to ovalbumin (CEELHAMLDVKPDADE)

Clones: 3ET-1C10

Format: Ascites fluid. Does not contain any preservative

Murine Isotype: IgG1, k

Applications:

- ELISA,
- Western blot
- Immunofluorescence

Recommended dilutions: 1/500-1/5000

Storage:

Store at -20°C (up to 3 years)
Avoid repeat freezing and thawing cycles. Does not contain any preservative

These products are intended for research purposes only.

MATERIAL SAFETY DATA SHEET INFORMATION

To the best of our knowledge, these products do not require a Material Safety Data Sheet. However, all the properties of these products (and, if applicable, each of their components) have not been thoroughly investigated. Therefore, we recommend that you use gloves and eye protection, and wear a laboratory coat when working with these products.