

## Product Data Sheet



<b>Product</b>	Chicken Whole Blood in Alsevers (1 : 2)
<b>Product No</b>	NB-43-0056-50
<b>Lot</b>	31-112-013KW27/23
<b>Date</b>	04/07/2023
<b>Volume</b>	50 ml (25 ml whole blood & 25 ml Alsevers)
<b>Presentation</b>	1 lot of 50 ml of chicken whole blood pooled from 2 female + 1 male animal
<b>Species of Origin</b>	Chicken
<b>Strain</b>	Crossbreed
<b>Sex</b>	female, male
<b>Remarks</b>	non medicated, non immunized
<b>Anticoagulant</b>	Alsever's solution
<b>Sterility</b>	aseptic
<b>Buffer</b>	Alsever's solution
<b>Preservative</b>	non
<b>Storage Conditions</b>	+2 - +8°C.
<b>Expiry Date</b>	18/07/2023
<b>Disclaimer</b>	<p>THIS MATERIAL SHOULD BE HANDLED AS IF CAPABLE OF TRANSMITTING INFECTIOUS AGENTS. PLEASE USE UNIVERSAL PRECAUTIONS.</p> <p>THIS PRODUCT IS BEING SOLD FOR RESEARCH AND, OR MANUFACTURING PURPOSES ONLY. IT IS NOT TO BE USED IN HUMANS OR ANIMALS AND FURTHER MANUFACTURING THAT WILL RESULT IN A FINISHED PRODUCT THAT CONTAINS VIABLE LEUKOCYTES IS PROHIBITED. FOR IN VITRO USE ONLY. THE USER ASSUMES ALL RESPONSIBILITY FOR ITS USAGE AND DISPOSAL, IN ACCORDANCE WITH ALL REGULATIONS.</p>
<b>Country of Origin</b>	Germany