



Normal Guinea pig serum blocking buffer

10% Normal Guinea Pig Serum

NB-23-00186-1 (100 ml)

NB-23-00186-2 (18 ml)



Normal Guinea pig serum blocking buffer

NB-23-00186-1; NB-23-00186-2

Storage: 2-8°C

INTENDED USE:

Normal Guinea Pig Serum Blocking Buffer is used in immunohistochemical staining procedures to block nonspecific staining when the host species of the secondary antibody is guinea pig.

KIT COMPONENTS:

Catalog No.	Size
NB-23-00186-2	18 mL (Ready to use)
NB-23-00186-1	100 mL (Ready to use)

Solution contains 0.05% NaN₃.

RECOMMENDED PROTOCOL:

1. Apply 100uL or enough volume of **Normal Guinea Pig Serum Blocking Buffer** on tissue before adding primary antibody and incubate for 10 minutes at room temperature.
2. Blot off Normal Guinea Pig Serum Blocking Buffer Solution (Do not rinse) and add primary antibody.

PRECAUTION:

Please wear gloves and take other necessary precautions.

FOR RESEARCH USE ONLY

Neo-Biotech
183, avenue Georges Clémenceau – 92000 Nanterre