

Product datasheet MON4103

MONOSAN[®]

Mouse anti-CD81, clone 1D6 (Monoclonal)

Clone no. 1D6

MONOSAN

Product name	Mouse anti-CD81, clone 1D6 (Monoclonal)
Host	Mouse
Applications	WB,IHC-P,IP,FC,IF
Species reactivity	human,sheep,chimpanzee,goat
Conjugate	Purified
Immunogen	OCI-LY8 cells aggregated by 5A6 (another CD81 antibody)
Isotype	IgG1
Clonality	Monoclonal
Clone number	1D6
Size	0.2 mg
Concentration	1.0 mg/ml
Format	-
Storage buffer	PBS with 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD81, clone 1D6 (Monoclonal)

Clone no. 1D6

MONOSAN

Additional info

Mouse anti Human CD81 antibody, clone 1D6 recognizes human CD81, a 26 kDa cell surface antigen also known as TAPA-1, and a member of the tetraspanin family. CD81 is widely expressed on human leucocytes and appears to be involved in a variety of cellular leucocytes including activation, proliferation and differentiation. Mouse anti Human CD81 antibody, clone 1D6 is a potent CD81 reagent, induces homotypic adhesion and has powerful anti-proliferative effects.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES