

NeoFuge M12 - Table Top Genius Micro-Centrifuge (Programmable)

#Cat: NB-12-4002 Size: 1unit

Conceptualized and developed to provide optimum performance along with exceptional ease-of-use, NeoFuge M12 begins where others end. Its state-of-the-art construction optimizes performance as well as durability by offering a long maintenance-free long life. With unmatched features designed to provide best-in-class protection from accidents, NeoFuge M12 doesn't skimp when it comes to safety also. Change the way your work with the brilliance of NeoFuge M12.



Product Features:







Dual Shell Design And Tested Explosion Proof Construction

> Safer Low Voltage Device (24V) With Universal Power Supply (with Range Of 95 – 230 V 50/60 Hz)





Lid Lock Safety Feature Opens The Lid Automatically On Run Completion

> Automatic Internal Diagnosis And Error Display On Screen For Unmatched Convenience.





Intelligent Imbalance Detection Feature Automatically Switches Off The Device In Case Of Imbalance While The Device Is Running.

12 X 1.5/2ml closed rotor





Microprocessor Controlled For Greater Efficiency





Variable Speed Setting From 500 RPM To 15000 RPM With Increments Of 100 RPM Or As Per Desired RCF.

> Intuitive Simple Interface For Quick And Convenient Setting With Last Run Memory Feature.





Variable Timer Setting From 1 To 999 Minutes With Min-sec Display And Continuous Mode Function.

> Unprecedented Convenience With One Touch Short Spin And One Touch Toggle Between RCF & RPM.





Short spin feature for short centrifugation & fast pelleting



Specifications	NeoFuge M12
Rotor capacity	12 x 1.5 ml / 2.0 ml micro tubes
Maximum RPM / RCF	15000 / 15596 g
Display	Large back-it LCD display
Imbalance detection	Yes
Noise level	< 60 db
Rotor	Rotor with Screw on lid (for better contamination
	protection)
Motor type	BLDC – Brushless DC Motor
Programmable	Yes, USB connectivity for data log & remote operation
	programming
Timer setting	Countdown timer from 1-999 mins (min-sec mode) &
	infinite
	mode
Dimension (WxDxH)	262 x 230 x 131 mm
Rotors and adaptors (included in	12 place x 1.5/2.0 ml
standard pack)	Reduction adaptors for 0.4/0.5 ml microtubes
Optional rotors	PCR strip rotor (2x8x0.2 ml with metal safety lid)

Ordering information:

Cat Number	Details
NB-12-4002	NeoFuge M12