

Multi Microplate Shaker #Cat: NB-12-0017

This shaker is a high performance microplate shaker which accommodates up to four microplates. Shaking speed and time are programmable via the keypad, while the status parameters are displayed on the LED in real time via an easy to use interface.

Features:

- Variable speed between 200 to 1,500 rpm
- DC brushless motor for gentle and powerful motion
- Orbital of 3mm for efficient mixing in small sample vessels
- Continuous or timed operation
- Simple spring holders for easy and speed of plate installation



Specifications:

Cat#	NB-12-0017
Speed	300-1400 rpm/min Standard
	300-1000 rpm/min Deepwell
Timer	59h99min
Orbit	3 mm
Weight (kg)	7.0
Power	100~240V 40W
Dimensions	280 × 270 × 110 mm