

PPP1R9A antibody

Product Information

Catalog No.:	FNab06710
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds to actin filaments(F-actin) and shows cross-linking activity. Binds along the sides of the F-actin. May be involved in neurite formation. Inhibits protein phosphatase 1-alpha activity(By similarity).

Immunogen information

Immunogen:	protein phosphatase 1, regulatory(inhibitor) subunit 9A
Synonyms:	Neurabin-1 Neurabin-I Neural tissue-specific F-actin-binding protein I Protein phosphatase 1 regulatory subunit 9A PPP1R9A KIAA1222
Observed MW:	180 kDa
Uniprot ID :	Q9ULJ8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab06710(PPP1R9A antibody) at dilution of 1:300