




- **Autofil[®] 2** is available as a complete system or funnel only
- Funnel is compatible with GL45 thread bottles
- **Autofil[®] 2** units include adapter for connecting to vacuum hose
- **Autofil[®] 2** design includes the **NEW EZlabpure[™]** ergonomic bottle shape with molded finger grip
- **Autofil[®] 2** patented GL45 caps provides a convenient gripping surface with a frosted area for labeling
- Made with USP Class VI materials
- Manufactured in a Class 100,000 clean-room
- Manufactured in an ISO 13485 Facility in North America
- Optional, patented docking hardware sold separately



The patented, next generation **Autofil[®] 2** features the **NEW High Flow Foxx Velocity 0.1µm, 0.2µm, and 0.45µm PES membranes**. Perfect for filtering FBS and other viscous solutions. Our membrane is low protein binding and surfactant-free. **Autofil[®] 2** is available in 250mL, 500mL, and 1000mL complete units with the **NEW EZlabpure[™]** media bottle or funnel only bottle-top filters. All units are sterile and manufactured from USP Class VI polystyrene in an ISO 13485 clean room, FDA regulated facility.


Autofil[®] 2 is compatible with the patented **Autofil[®] Pedestal** and **NEW Cradle Ring 2**. These optional docking stations allow for hands free filtering and provide a stationary vacuum hose connection, allowing for efficient set-up and quick use while preventing spills and tissue culture hood contamination.

- 1 Patented single-piece design made entirely from polystyrene material
- 2 Features **NEW High Flow Foxx Velocity 0.2µm PES membrane**
- 3 Forward-facing, large, raised graduation marks for easy reading
- 4 Easy to read label makes identifying membrane type and size easy when removing from packaging
- 5 **EZlabpure[™]** PS sterile media bottle with GL45 thread type



ISO 13485

Distributed by:
CliniSciences Group



	Item Number	Volume	Membrane Pore Size	Filter Diameter	Description
Full Assembly	116-71132-FLS	250mL	0.1µm PES	66mm	250mL Full Assembly 0.1µm PES, 12/pk
	116-81132-FLS	500mL	0.1µm PES	90mm	500mL Full Assembly 0.1µm PES, 12/pk
	116-91132-FLS	1000mL	0.1µm PES	90mm	1000mL Full Assembly 0.1µm PES, 12/pk
	11012-FLS	250mL	0.2µm PES	66mm	250mL Full Assembly 0.2µm PES, 12/pk
	11022-FLS	500mL	0.2µm PES	90mm	500mL Full Assembly 0.2µm PES, 12/pk
	11032-FLS	1000mL	0.2µm PES	90mm	1000mL Full Assembly 0.2µm PES, 12/pk
	11412-FLS	250mL	0.45µm PES	66mm	250mL Full Assembly 0.45µm PES, 12/pk
	11422-FLS	500mL	0.45µm PES	90mm	500mL Full Assembly 0.45µm PES, 12/pk
	116-3233-RLS	1000mL	0.45µm PES	90mm	1000mL Full Assembly 0.45µm PES, 12/pk
Funnel Only	116-11132-FLS	250mL	0.1µm PES	66mm	250mL Funnel Only 0.1µm PES, 12/pk
	116-21132-FLS	500mL	0.1µm PES	90mm	500mL Funnel Only 0.1µm PES, 12/pk
	116-31132-FLS	1000mL	0.1µm PES	90mm	1000mL Funnel Only 0.1µm PES, 12/pk
	11512-FLS	250mL	0.2µm PES	66mm	250mL Funnel Only 0.2µm PES, 12/pk
	11522-FLS	500mL	0.2µm PES	90mm	500mL Funnel Only 0.2µm PES, 12/pk
	11532-FLS	1000mL	0.2µm PES	90mm	1000mL Funnel Only 0.2µm PES, 12/pk
	11612-FLS	250mL	0.45µm PES	66mm	250mL Funnel Only 0.45µm PES, 12/pk
	11622-FLS	500mL	0.45µm PES	90mm	500mL Funnel Only 0.45µm PES, 12/pk
	11632-FLS	1000mL	0.45µm PES	90mm	1000mL Funnel Only 0.45µm PES, 12/pk
Optional Accessories	11712-FLS	250mL	-	-	250mL EZlabpure™ Media Bottle, PS, 12/pk
	11722-FLS	500mL	-	-	500mL EZlabpure™ Media Bottle, PS, 12/pk
	11732-FLS	1000mL	-	-	1000mL EZlabpure™ Media Bottle, PS, 12/pk
	1124-FLS	-	-	-	Autofil® Quick-Connect Cradle Ring 2, 1/pk
	1116-RLS	-	-	-	Autofil® Quick-Connect Pedestal, 1/pk



Autofil® is a registered trademark of Foxx Life Sciences.

CliniSciences Group

Austria



Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>

Belgium



Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>

Denmark



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: danmark@clinisciences.com
Web: <https://www.clinisciences.com>

Finland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>

France



Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>

Germany



Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mausier-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>

Iceland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>

Ireland



Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>

Italy



Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>

Netherlands



Company: CliniSciences B.V.
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>

Norway



Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>

Poland



Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>

Portugal



Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>

Spain



Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>

Sweden



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>

Switzerland



Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>

UK



Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
IWeb: <https://www.clinisciences.com>

USA



Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>