



THE ELISA EXPERTS



Distributed by:
CliniSciences Group



Distributed by:
CliniSciences Group

THE ELISA EXPERTS

Company profile

CANDOR Bioscience GmbH is an internationally oriented, innovative, profitable and privately owned enterprise. CANDOR develops, produces and distributes a comprehensive selection of premium solutions for immunoassays.

CANDOR offers more than 50 different products including optimizers, blockers, stabilizers and buffer solutions for immunoassays. By using these solutions the reliability of results can be

improved, the assay can be simplified and process times can be reduced.

More than 600 customers in 35 countries use CANDOR products and benefit from the CANDOR expertise in assay development and sample measurement. The CANDOR experts are well experienced in research, development and production and can interact with laboratory managers at their own level.

CANDOR offers its customers high quality across the board.

CANDOR quality comprises highest product quality, highest quality standards and demanding quality control in production in addition to comprehensive product support and customer service. This exceptional quality standard can be guaranteed by our well-educated, experienced and motivated staff members.



Since the inception of CANDOR in 2004, production volume and sales have grown dynamically. From 2005 to 2010 the production capacity increased more than 120 fold. Each year has seen double-digit sales growth and the turnover in 2010 was 25 fold higher than in 2005. We expect dynamic growth to continue for years to come.

Distributed by:
CliniSciences Group

PERFECTION

Production & Customer Support

All CANDOR products are manufactured in the new manufacturing plant in Wangen, Germany. The products are developed and produced in our facility equipped with modern and high quality laboratory and production technology. When it comes to quality we always go one step further than standard quality specifications in our industry to enable the best product quality possible and a premium lot-to-lot consistency. CANDOR also runs an up-to-date industrial water-purification system for ultra-pure water in the new production facility with a capacity of 600 litres per hour.

CANDOR is interested in long-term business relationships and we are always committed to our customers' needs. The comprehensive product range meets every need for premium solutions for immunoassays. Furthermore, CANDOR offers individual product support along with fast, flexible and reliable order processing. Users receive detailed product information for application of CANDOR products.



Our intensive service also includes reliable, informed and individual care of customer requests concerning immunoassays. Due to our extensive ELISA experience we will always be available to help you with your challenges or questions.

Our customers trust in CANDOR products and the extraordinary CANDOR quality over the long term.



Distributed by:
CliniSciences Group

OPTIMIZERS, BLOCKERS, STABILIZERS

The CANDOR portfolio

Results of immunoassays should always be correct and reliable. Interferences, matrix effects and degradation of standards or conjugates can occasionally result in bad results, termed as false-positives or false-negatives in diagnostics. False results can be a safety hazard for patients in human diagnostics and a risk for consumers in food diagnostics. Furthermore it can be an expensive problem in any research application. CANDOR products contribute significantly to this scientific challenge. False results can be significantly reduced by using CANDOR products. In many cases

our solutions also help to simplify the assay procedures and to reduce processing times. Several of our products offer multiple functions and advantages in one solution and are therefore much appreciated by users.

Our portfolio comprises optimizers, blockers, stabilizers and other high-quality buffer solutions for immunoassays. Among the premium CANDOR solutions, LowCross-Buffer acts as the leading product.



LowCross-Buffer is the first solution that combines several functionalities of diluents in one product. It is not only an assay diluent. Additionally it can reduce HAMA-interference and furthermore helps to avoid matrix effects. CANDOR calls this new class of advantageous and multifunctional solutions "Optimizers" for immunoassays. They are outstanding compared to other diluents and buffers. LowCross-Buffer is a kind of binding gate-keeper. The innovative buffer solution, used instead of common assay diluents, suppresses low to medium affinity binding in the assay and only allows high affinity binding reactions.



For challenging assays LowCross-Buffer is often the only diluent to effectively avoid interference.

Increasing need for reliability of immunological assays requires replacement of the usual diluents with LowCross-Buffer e.g. in drug development, detection of sensitive new biomarkers or in diagnostics to avoid false results. Thus there is a continuously growing demand of this innovative CANDOR buffer solution.

Worldwide more than 90 pharmaceutical companies and their CROs (Contract Research Organisations) regularly use LowCross-Buffer for assay development, validation and sample measurement. They can trust their results. LowCross-Buffer is a useful tool to reduce interference effects in immunoassays. This improves precision and accuracy and enables more reliable results. The use of LowCross-Buffer in the production of diagnostics is therefore increasing due to the avoidance of false positives and false negatives.

Distributed by:
CliniSciences Group

CONVENIENCE

One step incubation assays


LowCross-HRP is another product with several functionalities. This product enables convenient one step incubation assays without negatively affecting the reliability of the results. This is due to both the optimizing and stabilizing properties of the solution.



	1	2	3	4	5	6	7	8	9	10	11	12	Endpoint
A	0,001	0,001	0,001	0,001	0,001	0,001	0,001	0,001	0,001	0,001	0,001	0,001	0,001
B	0,891	0,911	0,898	0,898	0,898	0,898	0,898	0,898	0,898	0,898	0,898	0,898	0,911
C	0,892	0,893	0,901	0,901	0,901	0,901	0,901	0,901	0,901	0,901	0,901	0,901	0,911
D	0,901	0,902	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,902
E	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,902
F	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,902
G	0,897	0,899	0,001	0,001	0,897	0,899	0,899	0,899	0,899	0,899	0,899	0,899	0,902

DEVELOPMENT OF SALES

calculated as 100 % in 2005



In addition to the innovative optimizers, CANDOR also offers blockers based on modern chemistries and newly developed stabilizers. CANDOR is the technical leader in the area of immunoassay stabilizers. Antibody Stabilizer, Liquid Plate Sealer and HRP-Protector protect the native conformation of antibodies at room temperature or in the refrigerator much better than many competitor products.





Distributed by:
CliniSciences Group



EVERYTHING UNDER CONTROL

Production at CANDOR

All CANDOR products are manufactured in the new manufacturing plant in Wangen, Germany.

The products are developed, produced and controlled in our own facility equipped with modern and high quality laboratory and production technology.

The company is certified according to DIN EN ISO 9001:2008. All working steps are described in SOPs. The water purification system, the laboratory balances and the industrial scales meet highest international quality standards. Our company sets a high value on batch consistency and high product quality "Made in Germany". CANDOR always works to outperform the required quality standards to ensure the best reliability for users.



CANDOR products show the maximum shelf life at dispatch. The complete quality control ensures a constant product quality.

The premium packaging also contributes to this. The CANDOR quality is facilitated by our dependable team with well-educated, experienced and motivated staff members.

Distributed by:
CliniSciences Group

FAST AND FLEXIBLE

CANDOR Customer Service

CANDOR is interested in long-term business relationships and we are always committed to our customers' needs. All aspects of the company are strongly oriented to serve customers according to their specific needs.

CANDOR supplies very fast and flexible service. No one will be left on their own concerning requests or challenges around immunoassays. The CANDOR experts work successfully on new immunoassay reagents to further optimize assays and their results. The CANDOR experts can discuss with laboratory managers and assay developers from their point of view, as they were engaged with assay development for years for research use, pre-clinical/clinical studies as well as for diagnostics use. Our long time experience in the field of ELISA is available to all customers.



Many users trust in CANDOR products, CANDOR customer service and the extraordinary CANDOR quality. To ensure this outstanding quality we have a motivated and highly skilled team. Our business philosophy is to interact with customers and colleagues always in a fair, respectful and friendly way. Therefore all CANDOR staff members will deal with interested individuals in this way.



CliniSciences Group

Austria
Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>

Germany
Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>

Norway
Company: CliniSciences AS
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland
Company: CliniSciences AG
Address: Fracht Ost Flughafen
Kloten CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>

Belgium
Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Iceland
Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Poland
Company: CliniSciences sp.z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



UK
Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark
Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>

Ireland
Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>

Portugal
Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>

USA
Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>



Finland
Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>

Italy
Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>

Spain
Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



France
Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>

Netherlands
Company: CliniSciences B.V.
Address: Gaetano Martinolaan 85, 6229
CS Maastricht - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>

Sweden
Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>